

eMedCheck

User's guide for eMedCheck for Palm PDA

July 14, 2008

Institute for Systems Research

University of Maryland

Cooperative Agreement Number U50/CCU302718 from the CDC to NACCHO supported this publication. Its contents are solely the responsibility of the University of Maryland and the Advanced Practice Center for Public Health Emergency Preparedness and Response of Montgomery County, Maryland, and do not necessarily represent the official views of CDC or NACCHO.

© Copyright 2008 University of Maryland and Montgomery County APC. All rights reserved.

User's guide for eMedCheck

1) Loading the Program

- Go to **Applications** and click on **eMedCheck** (it must be present if it was installed properly, else follow the instructions in the “eMedCheck Installation Guide” to re-install the software).
- Now you may either choose to enter the data for a household (Go to Step 2) or verify the data saved in the database. To view the records, please hit the **View Records** button.

2) Input the number of people in the first family

- Click on the number in the numeric keypad below the question.

Note: i) *A maximum of 9 people is allowed to be entered per household. If there are more than 9 members, process the first 9 and the remaining in the next iteration.*

3) Enter Initials

- Enter the initials of the people in the household by tapping on the corresponding dotted lines.

Note: *Click on “abc” at the bottom of the keypad if an alphabetic keypad is not present on the display.*

- Press **Continue** when you are done entering the initials of all family members or press **Cancel** to change the number of people in the household.

4) Enter medical data

- Enter the details for all family members by checking the appropriate checkbox corresponding to each member.

Important Note:

- i) **Pregnant Or Under-8** -> Check this checkbox if the person is Pregnant or under the age of 8 or *if you're not sure of the answer.*
- ii) **Allergic to Doxy** -> Check this checkbox only if you're sure the person is allergic to Doxy (Vibramycin), Minocycline (Minocin) or Tetracycline (Achromycin, Brodspec, EmTet, Sumycin, Tetracap, Panmycin). If not sure, don't check this checkbox.

iii) Allergic to Cipro -> Check this checkbox only if you're sure the person is allergic to Ciprofloxacin (Cipro), Levofloxacin (Levaquin), Ofloxacin (Floxin), Gatifloxacin (Tequin) or Moxifloxacin (Avalox). If not sure, don't check this checkbox.

- Click **Continue** to continue or **Cancel** to start over.

5) Observe the results

- Observe the results displayed for each member of the family.

- Hit **Save** to save this information to the database.

- Else, hit **Restart** to start over without saving the results to the database.

6) Observe the database

- See that the latest information has been saved to the database.

- You may choose to either:

i) Continue to the next patient by clicking on **Next Patient**, Or

ii) Reset the database and set all counters to zero by clicking on **Reset Data**

Note: *You'll be prompted to confirm if you want to reset the database. Choose **Yes** if you're sure you want to erase all records.*